

Phospho-c-Jun-S63 Rabbit pAb

货号: B11136

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mouse embryonal fibroblast cell , 293T , Mus musculus , Homo sapiens , Gallus gallus , Chlorocebus aethiops , Rattus norvegicus , Oryza sativa , Sus scrofa
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	35kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic phosphorylated peptide around S63 of human c-Jun (NP_002219.1).
序列	LTSPD

靶点信息

研究背景	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
基因ID	3725

基因名	JUN
Swiss	P05412
别名	AP-1;AP1;c-Jun;jUN;c

产品验证



Western blot analysis of Phospho-c-Jun-S63 expressed in 293 using Phospho-c-Jun-S63 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn