

CTSK Rabbit pAb

货号: B11209

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Meretrix meretrix , Rattus norvegicus , Mus musculus , Homo sapiens IHC: Mus musculus IF: Mus musculus
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	36kDa
实测分子量	26KDa/37KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SK-MEL-28,U-87MG,MG-63
细胞定位	Lysosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 115-329 of human CTSK (NP_000387.1).
序列	APDSVDYRKKGYVTPVKNQGCQSCWAFSSVGALEGQLKKKTGKLLNLSPQNLVDCVSENDGCGGGYMTNAFQYVQK NRGIDSEDAYPYVGQEESCMYNPTGKAAKCRGYREIPEGNEKALKRAVARVGPVSVAIDASLTSFQFYSGVYYDESCNS DNLNHAVLAVGYGIQKGNKHWIIKNSWGENWGNKGYILMARNKNNACGIANLASFPKM

靶点信息

研究背景	The protein encoded by this gene is a lysosomal cysteine proteinase involved in bone remodeling and resorption. This protein, which is a member of the peptidase C1 protein family, is predominantly expressed in osteoclasts. However, the encoded protein is also expressed in a significant fraction of human breast cancers, where it could contribute to tumor invasiveness. Mutations in this gene are the cause of pycnodystosis, an autosomal recessive disease characterized by osteosclerosis and short stature.
基因ID	1513
基因名	CTSK
Swiss	P43235
别名	CTSK;CTS02;CTSO;CTSO1;CTSO2;PKND;PYCD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn