

Caspase-7 Rabbit pAb

货号: B11211

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Rattus norvegicus , Gallus gallus , gga , Mus musculus , Anatinae IF: Gallus gallus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	28kDa/31kDa/34kDa/37kDa
实测分子量	35KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human Caspase-7 (NP_001218.1).
序列	IYGKDGVTPIKDLTAHFRGDRCKTLLEKPKLFFIQACRGTELDG IQADSGPINDTDANPRYKIPVEADFLFAYSTVPGYYS WRSPGRGSWFVQALCSILE

靶点信息

研究背景	This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive pro enzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. The precursor of the encoded protein is cleaved by caspase 3 and 10, is activated upon cell death stimuli and induces apoptosis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	840
基因名	CASP7
Swiss	P55210
别名	CASP7;CASP-7;CMH-1;ICE-LAP3;LICE2;MCH3;caspase-7;Casp7

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn