

# EMR1 Rabbit pAb

货号: B11214

## 产品信息

|       |                                                                                                                  |
|-------|------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                  |
| 宿主    | Rabbit                                                                                                           |
| 克隆性   | Polyclonal                                                                                                       |
| 预测反应  | <b>WB:</b> Mus musculus<br><b>IHC:</b> Mus musculus , mouse<br><b>IF:</b> Rattus norvegicus , Mus musculus , mmu |
| 应用    | WB IF/ICC                                                                                                        |
| 推荐浓度  | <b>WB:</b> 1:100 - 1:500<br><b>IF/ICC:</b> 1:50 - 1:200                                                          |
| 理论分子量 | 78kDa/82kDa/90kDa/95kDa/97kDa                                                                                    |
| 实测分子量 | 160KDa                                                                                                           |
| 形式    | Liquid                                                                                                           |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.             |
| 偶联物   | Unconjugated                                                                                                     |
| 阳性对照  | Mouse lung                                                                                                       |
| 细胞定位  | Cell membrane,Multi-pass membrane protein                                                                        |
| 纯化    | Affinity purification                                                                                            |


## 抗原信息

|      |                                                                                                                             |
|------|-----------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 477-596 of human EMR 1 (NP_001965.3).         |
| 序列   | GVAFVSFVGMESVLNERFFKDHQAPLTTSEIKLMNSRVVGGIMTGEKKDGFSDPIIYTLENIQPKQKFERPICVSWSTDVK<br>GGRWTSFGCVILEASETYTICSCNQMANLAVIMASGEL |

## 靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | This gene encodes a protein that has a domain resembling seven transmembrane G protein-coupled hormone receptors (7TM receptors) at its C-terminus. The N-terminus of the encoded protein has six EGF-like modules, separated from the transmembrane segments by a serine/threonine-rich domain, a feature reminiscent of mucin-like, single-span, integral membrane glycoproteins with adhesive properties. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. |
| 基因ID  | 2015                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 基因名   | ADGRE1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Swiss | Q14246                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| 别名    | ADGRE1;EMR1;TM7LN3;F4/80                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)