

SLC5A1 Rabbit pAb

货号: B11239

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Sus scrofa , Rattus norvegicus IHC: Gallus gallus
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	60kDa/73kDa
实测分子量	73KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human SLC5A1 (NP_000334.1).
序列	ECEKYCGTKVGCTNIAYPTLVVELMPNGLRGLMLSVMLASLMSSLTSIFNSASTLFTMDIYAKVRKRASEKELMIAGRLFILVLIGISIAWVPIVQSAQSG

靶点信息

研究背景	This gene encodes a member of the sodium-dependent glucose transporter (SGLT) family. The encoded in tegral membrane protein is the primary mediator of dietary glucose and galactose uptake from the intestinal lumen. Mutations in this gene have been associated with glucose-galactose malabsorption. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	6523

基因名	SLC5A1
Swiss	P13866
别名	SLC5A1;D22S675;NAGT;SGLT1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn