

MYH4 Rabbit pAb

货号: B11322

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Sus scrofa , Mus musculus IF: Rattus norvegicus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	250KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse skeletal muscle,Rat skeletal muscle
细胞定位	Cytoplasm,myofibril
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1450-1550 of human MYH4 (NP_060003.2).
序列	QRNFDKVLAEWKQKYEETQAELEASQKESRSLSTELFKVKKNAYEESLDHLETLKRENKNLQQEISDLTEQIAEGGKHIHELEKVKKQLDHEKSELQTSLEE

靶点信息

研究背景	Muscle contraction.
基因ID	4622

基因名	MYH4
Swiss	Q9Y623
别名	MYH4;MYH2B;MyHC-2B;MyHC-IIb;myosin-4

产品验证



Western blot analysis of MYH4 expressed in Mouse skeletal muscle,Rat skeletal muscle using MYH4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn