

ZO-1 Rabbit pAb

货号: B11356

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Pig cells , Mus musculus , Sus scrofa IF: Gallus gallus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	186kDa/195kDa
实测分子量	250KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse heart
细胞定位	Cell junction,Cell membrane,Cytoplasmic side,Peripheral membrane protein,gap junction,tight junction
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 18-299 of human ZO-1 (NP_003248.3).
序列	IWEQHTVTLHRAPGFGFGIAISGGRDNPHFQSGETSIVISDVLKGGPAEGQLQENDRVAMVNGVSMDNVEHAFVQQLR KSGKNAKITIRKKKVQIPVSRPDPEPVS DNEEDSYDEEIHDPGRSGVGNRRSEKIWPRDRSASRERSLSPRSDRRSVA SSQPAKPTKVTLVKSRKNEEYGLRLASHIFVKEISQDSLAARDGNIQEGDVVLKINGTVTENMSLTDAKTLIERSKGKLMV VQRDERATLLNVPDLSDSIHSANASERDDISEIQSLASDH

靶点信息

研究背景	This gene encodes a protein located on a cytoplasmic membrane surface of intercellular tight junctions. The encoded protein may be involved in signal transduction at cell-cell junctions. Alternative splicing of this gene results in multiple transcript variants.
------	---

基因ID	7082
基因名	TJP1
Swiss	Q07157
别名	TJP1;ZO-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn