

TGF beta Receptor I (TGFBR1) Rabbit pAb

货号: **B11399**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Rattus norvegicus , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	56kDa
实测分子量	56KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse lung
细胞定位	Cell junction,Cell membrane,Single-pass type I membrane protein,tight junction
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TGF beta Receptor I (TGFBR1) (NP_004603.1).
序列	CCNQDHCNKIELPTTVKSSPGLGPVELAAVIAGPVCFCISLMLMVYICHNRTVIHHRVPNEEDPSLDRPFISEGTTLKDLIYD MTTSGSGSGLPLLVRT

靶点信息

研究背景	The protein encoded by this gene forms a heteromeric complex with type II TGF-beta receptors when bound to TGF-beta, transducing the TGF-beta signal from the cell surface to the cytoplasm. The encoded protein is a serine/threonine protein kinase. Mutations in this gene have been associated with Loeys-Dietz aortic aneurysm syndrome (LDAS). Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	7046

基因名	TGFBR1
Swiss	P36897
别名	TGFBR1;AAT5;ACVRLK4;ALK-5;ALK5;ESS1;LDS1;LDS1A;LDS2A;MSSE;SKR4;TGFR-1;tbetaR-I;TBRI;TBR-i

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn