

LAMP1 Rabbit pAb

货号: B11400

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens IF: Homo sapiens , Rattus norvegicus , Mus musculus IHC: Mus musculus
应用	WB IF/ICC IHC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200 IHC: 1:50 - 1:200
理论分子量	38kDa/44kDa
实测分子量	45KDa/110KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,Mouse kidney,Rat brain
细胞定位	Cell membrane,Endosome membrane,Late endosome,Lysosome membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human LAMP1 (NP_005552.3).
序列	DKKYRCVSGTQVHMNNVTVTLHDATIQAYLSNSSFSRGETRCEQDRPSPTTAPPAPPSPSPSPVPKSPSVDKYNVSGTNGTCLLASMGQLNLTYERKDNT

靶点信息

研究背景	The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.
------	---

基因ID	3916
基因名	LAMP1
Swiss	P11279
别名	LAMP1;CD107a;LAMP A;LGP120

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn