

Phospho-RIP3-T231/S232 Rabbit pAb

货号: B11414

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus ICH: Mus musculus IF: Mus musculus ELISA: Mus musculus
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	53kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	L929
细胞定位	cytoplasm,cytosol,nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic phosphorylated peptide around T231 & S232 of mouse Ripk3 (NP_064339.2).
序列	DKTSLI

靶点信息

研究背景	The product of this gene is a member of the receptor-interacting protein (RIP) family of serine/threonine p rotein kinases, and contains a C-terminal domain unique from other RIP family members. The encoded pr otein is predominantly localized to the cytoplasm, and can undergo nucleocytoplasmic shuttling depende nt on novel nuclear localization and export signals. It is a component of the tumor necrosis factor (TNF) r eceptor-I signaling complex, and can induce apoptosis and weakly activate the NF-kappaB transcription fa ctor.
------	--

基因ID	56532
基因名	
Swiss	Q9QZL0
别名	R; Rip3; AW107945; 2610528K09Rik

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn