

Hsp90α Rabbit pAb

货号: **B11417**

产品信息

反应	Human,Mouse,Rat,Monkey
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Mouse Embryonic Stem Cells WB: Human liver cell , 293T , Homo sapiens , Mus musculus IF: Human liver cell , Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:20 - 1:50
理论分子量	84kDa/98kDa
实测分子量	90KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,NIH/3T3,Rat brain,PC-12
细胞定位	Cell membrane,Cytoplasm,Melanosome
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 629-732 of human Hsp90α (NP_005339.3).
序列	AAKKHLEINPDHSIIETLRQKAEADKNDKSVKDLVILLYETALLSSGFSLEDPQTHANRIYRMIKGLGLIDEDDPTADDTSAAVTEEMPPLGDDDDTSRMEEVD

靶点信息

研究背景	The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	3320
基因名	HSP90AA1
Swiss	P07900
别名	HSP90AA1;EL52;HEL-S-65p;HSP86;HSP89A;HSP90A;HSP90N;HSPC1;HSPCA;HSPCAL1;HSPCAL4;HSPN;Hsp103;Hsp89;Hsp90;LAP-2;LAP2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn