

THBS2 Rabbit pAb

货号: **B11438**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens , Mus musculus WB: Homo sapiens , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	129kDa
实测分子量	150kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	basement membrane,extracellular region
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 20-280 of human THBS2 (NP_003238.2).
序列	HQDKDTTFDLFSISNINRKTIGAKQFRGPDGPVPAYRFVRFDYIPPVNADDLSKITKIMRQKEGFFLTAQLKQDGKSRGTLL ALEGPGLSQRQFEIVSNGPADTLDLTYWIDGTRHVVSLLEDVGLADSQWKNVTVQVAGETYSLHVGCDLIDSFALDEPFYE HLQAEKSRMYVAKGSARESHFRGLLQNVHLVFENSVEDILSKKGCQQGQGAEINAISENTETLRLGPHVTTEYVGPSSERR PEVCERSCEELGNMVQEL

靶点信息

研究背景	The protein encoded by this gene belongs to the thrombospondin family. It is a disulfide-linked homotrimeric glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. This protein has been shown to function as a potent inhibitor of tumor growth and angiogenesis. Studies of the mouse counterpart suggest that this protein may modulate the cell surface properties of mesenchymal cells and be involved in cell adhesion and migration.
基因ID	7058
基因名	THBS2
Swiss	P35442
别名	THBS2;TSP2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn