

TET2 Rabbit pAb

货号: B11448

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Mus musculus , Sus scrofa , Homo sapiens WB: Mus musculus , human IHC: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	130kDa/133kDa/223kDa
实测分子量	280KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 303-470 of human TET2 (NP_001120680.1).
序列	NASKLAAMLNTCSFQKPEQLQQQKSVFEICPSPAENNIQGTTKLASGEEFCSGSSSNLQAPGGSSERYLKQNEMNGAYFK QSSVFTKDSFSATTTTPPPSQLLLSPPPPLPQVPQLPSEGKSTLNGGVLEEHHHYPNQSNTTLLREVKIEGKPEAPPSQSPN PSTHVC

靶点信息

研究背景	The protein encoded by this gene is a methylcytosine dioxygenase that catalyzes the conversion of methylcytosine to 5-hydroxymethylcytosine. The encoded protein is involved in myelopoiesis, and defects in this gene have been associated with several myeloproliferative disorders. Two variants encoding different isoforms have been found for this gene.
基因ID	54790
基因名	TET2
Swiss	Q6N021
别名	MDS;KIAA1546;TET2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn