

Hsp90β Rabbit pAb

货号: B11506

产品信息

反应	Human,Mouse,Rat,Monkey
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Human liver cell WB: Human liver cell , Rat cells , Mus musculus , Homo sapiens , Rattus norvegicus , Gallus gallus IHC: Mouse Embryonic Stem Cells
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	83kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,DU145,SW620,COS-7,NIH/3T3,HeLa
细胞定位	Cytoplasm,Melanosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Hsp90β (NP_031381.2).
序列	MPEEVHHGEEEVETFAFQAEIAQLMSLIINTFYSNKEIFLRELISNASDALDKIRYESLTDP SKLDSGKELKIDIIPNPQERTLTLVDTGIGMTKADLINNLGTIAKSGTKAFMEALQAGADISMIGQFGVGFYSAYLVAEKVVVITKHNDDEQYAWESSAGGSFTVRADHGEPIGRGTKVILHLKEDQTEYLEERRVKEVVKKHSQFIGYPITLYLEKEREKEISDDEAEEEEKGEKEEEDKDDEEKPKIEDVGSDEEDDSGKDKKKKTKKIKEYIDQEELNKTPIWTRNPDDITQEE

靶点信息

研究背景	This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes.
基因ID	3326
基因名	HSP90AB1
Swiss	P08238
别名	D6S182;HSP84;HSP90B;HSPC2;HSPCB;Hsp90 beta;HSP90AB1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn