

KPNA3 Rabbit pAb

货号: B11513

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Human hepatocellular carcinoma cell line , Homo sapiens , Mus musculus , Sus scrofa Co-IP: Sus scrofa
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	57kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,Jurkat,A-549,U-251MG,HeLa
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human KPNA3 (NP_002258.2).
序列	MAENPSLENHRIKSFKNKGRDVTMRRHRNEVTVELRKNKRDEHLLKKNVNPQEESESDVDADFKAQNVTLQNA TSDNPVVQLSAVQAARKLLSSDRNPPIDDLIKSGILPILVKCLERDDNPSLQFEAAWALTNIASGTSAGTQAVVQSNVPLF LRLLRSPHQNVCEQAVWALGNIIIGDGPQCRDYGVISLGVVKPLLSFISP

靶点信息

研究背景	The transport of molecules between the nucleus and the cytoplasm in eukaryotic cells is mediated by the nuclear pore complex (NPC), which consists of 60-100 proteins. Small molecules (up to 70 kD) can pass through the nuclear pore by nonselective diffusion while larger molecules are transported by an active process. The protein encoded by this gene belongs to the importin alpha family, and is involved in nuclear protein import.
基因ID	3839
基因名	KPNA3
Swiss	O00505
别名	KPNA3;IPOA4;SRP1;SRP1gamma;SRP4;hSRP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn