

ZEB2 Rabbit pAb

货号: B11530

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus IHC: Mus musculus
应用	WB IHC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200
理论分子量	133kDa/136kDa
实测分子量	210KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1050-1214 of human ZEB2 (NP_055610.1).
序列	SGEKPYQCDKCGKRFSHSGSYSQHMHRYSYCKREAEEREAAEREAREKGHLEPTTELLMNRAYLQSITPQGYSDSEERESMPRDGESEKEHEKEGEDGYGKLGRQDGDDEEEEEEESENKSMDDTDPETIRDEEETGDHSMDDSSDGGKMETKSDHEEDNMEDGM

靶点信息

研究背景	The protein encoded by this gene is a member of the Zfh1 family of 2-handed zinc finger/homeodomain proteins. It is located in the nucleus and functions as a DNA-binding transcriptional repressor that interacts with activated SMADs. Mutations in this gene are associated with Hirschsprung disease/Mowat-Wilson syndrome. Alternatively spliced transcript variants have been found for this gene.
------	--

基因ID	9839
基因名	ZEB2
Swiss	O60315
别名	ZEB2;HSPC082;SIP-1;SIP1;SMADIP1;ZFHX1B

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn