

SLC32A1 Rabbit pAb

货号: B11552

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	57kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spinal cord,Mouse brain
细胞定位	Cytoplasmic vesicle membrane,Multi-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-525 of human SLC32A 1 (Q9H598).
序列	MATLLRSKLSNVATSVSNKSQAKMSGMFARMGFGAATDEEAVGFAHCDDLDFEHRQGLQMDILKAEGEPCGDEGAEA PVEGDIHYQRGSGAPLPPSGSKDQVGGGGEGFGGHDKPKITAWAGWNVTNAIQGMFVLGLPYAILHGGYLGLFLIIFAAV VCCYTGKILIACLYEENEDGEVVRVSDSYVAIANACCAPRFPTLGGRVVNVAQIIELVMTCILYVVVSGNLMYNSFPGLPVSQ KWSIIATAVLLPCAFLKNLKAVSKFSLLCTLAHFVINILVIAYCLSRARDWAWEKVKFYIDVKKFPISIGIIVFSYTSQIFLPSLE GNMQQPSEFHCMMNWITHIAACVLKGLFALVAYLTWADETKEVITDNLPGSIRAVVNIFLVAKALLSYPLPFFAAVEVLEKSL FQEGSRAFFPACYSGDGRKSWGTLRLCALVVFTLLMAIYVPHFALLMGLTGSLTGAGLCFLLPSLFHLRLLWRKLLWHQV FFDVAIFVIGGICSVSGFVHSLEGLIEAYRTNAED

靶点信息

研究背景	The protein encoded by this gene is an integral membrane protein involved in gamma-aminobutyric acid (GABA) and glycine uptake into synaptic vesicles. The encoded protein is a member of amino acid/polyamine transporter family II.
基因ID	140679
基因名	SLC32A1
Swiss	Q9H598
别名	SLC32A1;VGAT;VIAAT

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn