

# Snail Rabbit pAb

货号: B11634

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>IHC:</b> Mouse Embryonic Stem Cells , Homo sapiens , Rattus norvegicus <b>IF:</b> Cervical cancer cells <b>WB:</b> Cervical cancer cells , Mouse kidney cortex cells , Homo sapiens
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:100
理论分子量	29kDa
实测分子量	29kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human Snai 1 (NP_005976.2).
序列	RTHTGEKPFSCPHCSRAFADRSNLR AHLQTHSDVKKYQCQACARTFSRMSLLHKHQESGCSGCPR


## 靶点信息

研究背景	The Drosophila embryonic protein snail is a zinc finger transcriptional repressor which downregulates the expression of ectodermal genes within the mesoderm. The nuclear protein encoded by this gene is structurally similar to the Drosophila snail protein, and is also thought to be critical for mesoderm formation in the developing embryo. At least two variants of a similar processed pseudogene have been found on chromosome 2.
------	--



基因ID	6615
基因名	SNAI1
Swiss	O95863
别名	SNAI1;SLUGH2;SNA;SNAH;SNAIL;SNAIL1;dJ710H13.1

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)