

Bmi1 Rabbit pAb

货号: B11656

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: NCI-N87 cells , Homo sapiens CHIP: Homo sapiens
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:50 - 1:100
理论分子量	36kDa
实测分子量	43KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,K-562,Mouse heart,Mouse thymus,Rat thymus,Rat testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Bmi1 (NP_005171.4).
序列	LKDYYTLMDIAYIYTWRRNGPLPLKYRVRPTCKRMKISHQRDGLTNAGELESDSGSDKANSPAGGIPSTSSCLPSPSTPVQSPHPQFPHISSTMNGTSNSP

靶点信息

研究背景	This gene encodes a ring finger protein that is major component of the polycomb group complex 1 (PRC1). This complex functions through chromatin remodeling as an essential epigenetic repressor of multiple r egulatory genes involved in embryonic development and self-renewal in somatic stem cells. This protein also plays a central role in DNA damage repair. This gene is an oncogene and aberrant expression is asso ciated with numerous cancers and is associated with resistance to certain chemotherapies. A pseudogen e of this gene is found on chromosome X. Read-through transcription also exists between this gene and t he upstream COMM domain containing 3 (COMMD3) gene.
基因ID	648
基因名	
Swiss	P35226
别名	BMI1;FLVI2/BMI1;PCGF4;RNF51;flvi-2/bmi-1;Bmi1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn