

FOSL1 Rabbit pAb

货号: B11665

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Rattus norvegicus , Mus musculus IHC: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	17kDa/29kDa
实测分子量	40KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,C6
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FOSL1 (NP_005429.1).
序列	MFRDFGEPGPSSGNGGGYGGPAQPAAQAQKQFHLVPSINTMSGSQELQWMVQPHFLGPSSYPRPLTYPQYSPPQP RPGVIRALGPPPGVRRRPCEQI

靶点信息

研究背景	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. Several transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	8061
基因名	FOSL1
Swiss	P15407
别名	FOSL1;FRA;FRA1;fra-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn