

ADRB3 Rabbit pAb

货号: **B11676**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus IHC: Homo sapiens IF: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	43kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,K-562
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ADRB3 (NP_000016.1).
序列	MAPWPHENSSLAPWPDLPNTANTSGLPGVPEAALAGALLALAVLATVGGNLLVIVAIWTPRLQTMNTNVFTSLAA ADLVMGLLVPPAATLALT

靶点信息

研究背景	The protein encoded by this gene belongs to the family of beta adrenergic receptors, which mediate catecholamine-induced activation of adenylate cyclase through the action of G proteins. This receptor is located mainly in the adipose tissue and is involved in the regulation of lipolysis and thermogenesis.
------	--

基因ID	155
基因名	ADRB3
Swiss	P13945
别名	ADRB3;BETA3AR

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn