

TEAD4 Rabbit pAb

货号: B11697

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens WB、CO-IP: Homo sapiens IHC: Homo sapiens
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	34kDa/44kDa/48kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,K-562,SW620,HeLa,Mouse lung,Mouse skeletal muscle,Mouse liver,Rat lung,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-250 of human TEAD 4 (NP_003204.2).
序列	AREIQAKLKDQAAKDKALQSMAAMSSAQIISATAFHSSMALARGPGRPAVSGFWQGALPGQAGTSHDVKPFSQQTYAV QPPLPLPGFESPAGPAPSPSAPPAPPWQGRSVASSKLWMLEFSAFLEQQQDPDTYNKHLFVHI

靶点信息

研究背景	This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is preferentially expressed in the skeletal muscle, and binds to the M-CAT regulatory element found in promoters of muscle-specific genes to direct their gene expression. Alternatively spliced transcripts encoding distinct isoforms, some of which are translated through the use of a non-AUG (UUG) initiation codon, have been described for this gene.
------	---

基因ID	7004
基因名	TEAD4
Swiss	Q15561
别名	TEAD4;EFTR-2;RTEF1;TCF13L1;TEF-3;TEF3;TEFR-1;hRTEF-1B

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn