

# CDC20 Rabbit pAb

货号: B11727

## 产品信息

|       |                                                                                                                   |
|-------|-------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                   |
| 宿主    | Rabbit                                                                                                            |
| 克隆性   | Polyclonal                                                                                                        |
| 预测反应  | <b>WB:</b> Homo sapiens , Mus musculus<br><b>IF:</b> Mus musculus,Homo sapiens                                    |
| 应用    | WB IHC IF/ICC IP                                                                                                  |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000<br><b>IHC:</b> 1:50 - 1:100<br><b>IF/ICC:</b> 1:50 - 1:100<br><b>IP:</b> 1:500 - 1:1000 |
| 理论分子量 | 54kDa                                                                                                             |
| 实测分子量 | 55KDa                                                                                                             |
| 形式    | Liquid                                                                                                            |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.              |
| 偶联物   | Unconjugated                                                                                                      |
| 阳性对照  | PC-3                                                                                                              |
| 细胞定位  | Cytoplasm,centrosome,cytoskeleton,microtubule organizing center,spindle pole                                      |
| 纯化    | Affinity purification                                                                                             |


## 抗原信息

|      |                                                                                                                                                                                                                                                                                          |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human CDC20 (NP_001246.2).                                                                                                                                                                        |
| 序列   | MAQFAFESDLHSLQLDAIPNAPPARWQRKAKEAAGPAPSPMRAANRSHSAGRTPGRTPGKSSSKVQTTPSKPGGDRY<br>IPHRSAAQMEVASFLLSKENQPENSQTPTKKEHQKAWALNLNGFDVEEAKILRLSGKPQNAPEGYQNRCLKVLYSQKATPG<br>SSRKTCRYIPSLPDRILDAPEIRNDYYLNLVDWSSGNVLAVALDNSVYLWSASSGDILQLLQMEQPGEYISSVAWIKENY<br>LAVGTSSAEVQLWDVQQQKRLRNMTSHSAR |

## 靶点信息

|       |                                                                                                                                                                                                                                                  |
|-------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes, nuclear movement prior to anaphase and chromosome separation. |
| 基因ID  | 991                                                                                                                                                                                                                                              |
| 基因名   | CDC20                                                                                                                                                                                                                                            |
| Swiss | Q12834                                                                                                                                                                                                                                           |
| 别名    | CDC20;CDC20A;bA276H19.3;p55CDC                                                                                                                                                                                                                   |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)