

Clusterin alpha chain Rabbit pAb

货号: B11734

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	32kDa/48kDa/52kDa/53kDa/57kDa
实测分子量	38kDa/57kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,U-251MG,Human plasma,Mouse plasma
细胞定位	Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle,Endoplasmic reticulum,Microsome,Mitochondrion membrane,Nucleus,Peripheral membrane protein,Secreted,chromaffin granule,cytosol,secretory vesicle
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 220-449 of human Clusterin alpha chain (NP_976084.1).
序列	FPKSRIVRSLMPFSPYEPLNFHAMFQPFLEMIHEAQAMDIIHFHSPAFQHPPTEFIREGDDDRTVCREIRHNSTGCLRMKDQCDKCREILSVDCSTNNPSQAKLRRELDLSLQVAERLTRYNELLSYQWKMLNTSSLLEQLNEQFNWVSRLANLTQGEDQYYLRVTTVASHTSDSDVPSGVTEVVVKLFDSDPITVTPVEVSRKNPKFMETVAEKALQEYRKKHREE

靶点信息

研究背景	The protein encoded by this gene is a secreted chaperone that can under some stress conditions also be found in the cell cytosol. It has been suggested to be involved in several basic biological events such as cell death, tumor progression, and neurodegenerative disorders. Alternate splicing results in both coding and non-coding variants.
基因ID	1191
基因名	CLU
Swiss	P10909
别名	CLU;AAG4;APO-J;APOJ;CLI;CLU1;CLU2;KUB1;NA1/NA2;SGP-2;SGP2;SP-40;TRPM-2;TRPM2;clusterin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn