

THBS1 Rabbit pAb

货号: B11830

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: unknow , Homo sapiens IHC: Cervus nippon
应用	WB IHC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200
理论分子量	120kDa/129kDa
实测分子量	160KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat lung
细胞定位	Endoplasmic reticulum,Sarcoplasmic reticulum
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 19-352 of human THBS1 (NP_003237.2).
序列	NRIPESGGDNSVFDIFELTGAARKGSGRRLVKGPDPPSPAFRIEDANLIPPVPDDKFQDLVDAVRAEKGFLLLASLRQMKK TRGTLALERKDHSGQVFSVVSNGKAGTLDLSLTVQGKQHVVSVVEEALLATGQWKSITLFVQEDRAQLYIDCEKMENAEL DVPIQSVFTRDLASIALRLIAKGGVNDNFQGVLQNVRFVFGTTPEDILRNKGCSSTSVLLTLDNNVVNGSSPAIRTNYIGH KTKDLQAICGISDELSSMVLELRGLRTIVTTLQDSIRKVTEENKELANELRRPPLCYHNGVQYRNNEEWTVDSCTECHCQN SVTICKKVS

靶点信息

研究背景	The protein encoded by this gene is a subunit of a disulfide-linked homotrimeric protein. This protein is an adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. This protein can bind to fibronectin, laminin, type V collagen and integrins alpha-V/beta-1. This protein has been shown to play roles in platelet aggregation, angiogenesis, and tumorigenesis.
基因ID	7057
基因名	THBS1
Swiss	P07996
别名	THBS1;THBS;THBS-1;TSP;TSP-1;TSP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn