

Bid Rabbit pAb

货号: B11843

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Human fibrosarcoma cells , Mus musculus , Bos taurus , Rattus norvegicus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	11kDa/14kDa/21kDa/26kDa
实测分子量	22KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Rat thymus
细胞定位	Cytoplasm,Mitochondrion membrane,Mitochondrion membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human Bid (NP_001187.1).
序列	MDCEVNNGSSLRDECITNLLVFGFLQSCSDNSFRRELDALGHELPLAPQWEGYDELQTDGNRSSHSRLGRIEADSESQ EDIIRNIARHLAQVGDSMDRSIPPGLVNGLALQLRNTSRSEEDRNRDLATALEQLLQAYPRDMEKEKTMLVLALLLAKKVAS HTPSLLRDVFHTTVNFINQNLRTYVRSLARNGMD

靶点信息

研究背景	This gene encodes a death agonist that heterodimerizes with either agonist BAX or antagonist BCL2. The encoded protein is a member of the BCL-2 family of cell death regulators. It is a mediator of mitochondria l damage induced by caspase-8 (CASP8); CASP8 cleaves this encoded protein, and the COOH-terminal par t translocates to mitochondria where it triggers cytochrome c release. Multiple alternatively spliced trans cript variants have been found, but the full-length nature of some variants has not been defined.
------	--

基因ID	637
基因名	BID
Swiss	P55957
别名	BID;FP497

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn