

# F2R Rabbit pAb

货号: B11853

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>WB:</b> Danio rerio , Rattus norvegicus<br><b>IF:</b> Rattus norvegicus                           |
| 应用    | <a href="#">WB</a>                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000                                                                            |
| 理论分子量 | 47kDa                                                                                                |
| 实测分子量 | 70KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Mouse heart,Rat lung                                                                                 |
| 细胞定位  | Cell membrane,Multi-pass membrane protein                                                            |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 376-425 of human F2R ( NP_001983.2). |
| 序列   | SECQRYVYSILCKESSDPSSYNSSGQLMASKMDTCSSNLNNSIYKKLLT                                                                  |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Coagulation factor II receptor is a 7-transmembrane receptor involved in the regulation of thrombotic response. Proteolytic cleavage leads to the activation of the receptor. F2R is a G-protein coupled receptor family member. Alternative splicing results in multiple transcript variants. |
| 基因ID | 2149                                                                                                                                                                                                                                                                                           |

|       |                            |
|-------|----------------------------|
| 基因名   | F2R                        |
| Swiss | P25116                     |
| 别名    | F2R;CF2R;HTR;PAR-1;PAR1;TR |

产品验证



Western blot analysis of F2R expressed in Mouse heart,Rat lung using F2R Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)