

DLG1 Rabbit pAb

货号: B11882

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rattus norvegicus , Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	87-103kDa
实测分子量	114kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,LO2,HeLa,Jurkat
细胞定位	Basolateral cell membrane,Cell junction,Cell membrane,Endoplasmic reticulum membrane,Membrane,Peripheral membrane protein,postsynaptic cell membrane,postsynaptic density,sarcolemma,synapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human DLG1 (NP_004078.2).
序列	MPVRKQDTQRALHLL EYRSKLSQTEDRQLRSSIERVINIFQSNLFQALIDIQEFYEVTLTDNPKCIDRSKPSEPIQPVNTWEI SSLPSSTVTSETLPSSLSPSVEKYRYQDEDTPPQEHISPQITNEVIGPELVHVSEKNLSEIENVHGFVSHSHISPIKPT EAVLPS PPTVPVIPVLPVPAENTVILP

靶点信息

研究背景	This gene encodes a multi-domain scaffolding protein that is required for normal development. This protein may have a role in septate junction formation, signal transduction, cell proliferation, synaptogenesis and lymphocyte activation. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene, but the full-length nature of some of the variants is not known.
------	--

基因ID	1739
基因名	DLG1
Swiss	Q12959
别名	DLG1;DLGH1;SAP-97;SAP97;dj1061C18.1.1;hdlg

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn