

UVRAG Rabbit pAb

货号: B11885

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	35kDa/78kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa,Mouse liver,Rat liver,Rat brain
细胞定位	Chromosome,Early endosome,Endoplasmic reticulum,Late endosome,Lysosome,Midbody,centromere
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 501-699 of human UVRAG (NP_003360.2).
序列	KGHRKRASSENERLQYKTPPPSYNSALAQPVTTVPSMGETERKITSLSSSLDTSLDFSKENKKKGEDLVGSLNGGHANVH PSQEQGEALSGHRATVNGTLLPSEQAGSASVQLPGEFHPVSEALCCTVEQAEIIIGLEATGFASGDQLEAFNCIPVDSAV AVECDEQVLGEFEEFSRRIYALNENVSSFRRPRRSSDK

靶点信息

研究背景	This gene complements the ultraviolet sensitivity of xeroderma pigmentosum group C cells and encodes a protein with a C2 domain. The protein activates the Beclin1-PI(3)KC3 complex, promoting autophagy and suppressing the proliferation and tumorigenicity of human colon cancer cells. Chromosomal aberrations involving this gene are associated with left-right axis malformation and mutations in this gene have been associated with colon cancer.
------	--

基因ID	7405
基因名	UVRAG
Swiss	Q9P2Y5
别名	UVRAG;DHTX;VPS38;p63

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn