

Thioredoxin 1 (Trx1/TXN) Rabbit pAb

货号: B11893

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:500 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	9kDa/11kDa
实测分子量	12kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,H460,MCF7,HT-29,SKOV3
细胞定位	Cytoplasm,Nucleus,Secreted
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 20 to the C-terminus of human Thior edoxin 1 (Trx1/Thioredoxin 1 (Trx1/TXN)) (NP_003320.2).
序列	DKLVVVDFSATWCGPCKMIKPFHSLSEKYSNVIFLEVDVDDCQDVASECEVKCMPTFQFFKKGQKVGEFSGANKEKLE ATINELV

靶点信息

研究背景	The protein encoded by this gene acts as a homodimer and is involved in many redox reactions. The enc oded protein is active in the reversible S-nitrosylation of cysteines in certain proteins, which is part of the response to intracellular nitric oxide. This protein is found in the cytoplasm. Two transcript variants encod ing different isoforms have been found for this gene.
------	---

基因ID	7295
基因名	TXN
Swiss	P10599
别名	TXN;TRDX;TRX;TRX1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn