

# TEAD3 Rabbit pAb

货号: B11894

产品信息

|       |                                                                                                                                          |
|-------|------------------------------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                                          |
| 宿主    | Rabbit                                                                                                                                   |
| 克隆性   | Polyclonal                                                                                                                               |
| 预测反应  | <b>Co-IP:</b> Homo sapiens<br><b>WB:</b> Mus musculus<br><b>IP:</b> Mus musculus<br><b>IF:</b> Mus musculus<br><b>ChIP:</b> Mus musculus |
| 应用    | WB IHC IF/ICC                                                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200                                                     |
| 理论分子量 | 48kDa                                                                                                                                    |
| 实测分子量 | 49KDa                                                                                                                                    |
| 形式    | Liquid                                                                                                                                   |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.                                   |
| 偶联物   | Unconjugated                                                                                                                             |
| 阳性对照  | Mouse heart,Mouse lung,Rat heart                                                                                                         |
| 细胞定位  | Nucleus                                                                                                                                  |
| 纯化    | Affinity purification                                                                                                                    |


抗原信息

|      |                                                                                                                                                                                                                  |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 130-330 of human TEAD 3 (NP_003205.2).                                                                                             |
| 序列   | SSAQIVSASVLQNKFSPPSPLPQAVFSTSSRFWSSPPLLGGQQGPSQDIKPFAQPAYPIQPPLPPTLSSYEPLAPLPSAAASV<br>PVWQDRTIASSRLRLLEYSAFMEVQRDPDTYSKHLFVHIGQTNPAFSDPPLEAVDVRQIYDKFPEKKGGLKELYEGPPNA<br>FFLVKFWADLNSTIQEGPGAFYGVSSQYSSADSMTIS |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                              |
|-------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is predominantly expressed in the placenta and is involved in the transactivation of the chorionic somatomammotropin-B gene enhancer. Translation of this protein is initiated at a non-AUG (AUA) start codon. |
| 基因ID  | 7005                                                                                                                                                                                                                                                                                                                                                                         |
| 基因名   | TEAD3                                                                                                                                                                                                                                                                                                                                                                        |
| Swiss | Q99594                                                                                                                                                                                                                                                                                                                                                                       |
| 别名    | TEAD3;DTEF-1;ETFR-1;TEAD-3;TEAD5;TEF-5;TEF5                                                                                                                                                                                                                                                                                                                                  |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)