

FOXP3 Rabbit pAb

货号: B11931

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Rattus norvegicus
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	43kDa/44kDa/47kDa/49kDa
实测分子量	47KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse thymus
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 206-300 of human FOXP 3 (NP_054728.2).
序列	KVFEEPEDFLKHCQADHLLDEKGRAQCLLQREMVQSLEQQLVLEKEKLSAMQAHLAGKMALTKASSVASSDKGSCCIVA AGSQGPVVPAWSGPRE

靶点信息

研究背景	The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regul ators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively splice d transcript variants encoding different isoforms have been identified.
基因ID	50943

基因名	FOXP3
Swiss	Q9BZS1
别名	FOXP3;AIID;DIETER;IPEX;JM2;PIDX;XPID

产品验证



Western blot analysis of FOXP3 expressed in Mouse thymus using FOXP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn