

GSTK1 Rabbit pAb

货号: B11941

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens IP: Vitis amurensis Rupr WB: A. cantonensis
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:100
理论分子量	20kDa/24kDa/25kDa/31kDa
实测分子量	25KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,HeLa,HepG2,MCF7
细胞定位	Peroxisome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-226 of human GSTK1 (NP_057001.1).
序列	MGPLPRTVELFYDVLSPYSWLGFEILCRYQNIWNINLQLRP SLITGIMKDSGNKPPGLLPRKGLYMANDLKLLRHHLQIPIHF PKDFLSVMLEKGSLSAMRFLTAVNLEHPEMLEKASRELWMRVWSRNEDITEPQSILAAAEKAGMSAEQAQGLLEKIATPKV KNQLKETTEAACRYGAFGLPITVAHV DGGQTHMLFGSDRMELLAHLLGEKWMGPIPPAVNARL

靶点信息

研究背景	This gene encodes a member of the kappa class of the glutathione transferase superfamily of enzymes that function in cellular detoxification. The encoded protein is localized to the peroxisome and catalyzes the conjugation of glutathione to a wide range of hydrophobic substrates facilitating the removal of these compounds from cells. Alternative splicing results in multiple transcript variants.
基因ID	373156
基因名	GSTK1
Swiss	Q9Y2Q3
别名	GSTK1;GST;GST 13-13;GST13;GST13-13;GSTK1-1;hGSTK1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn