

ELN Rabbit pAb

货号: B11955

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	48kDa/56-69kDa
实测分子量	70KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Rat testis
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 27-200 of human ELN (NP_001265845.1).
序列	GGVPGAIPGGVPGGVFYPGAGLGALGGGALPGGKPLKPVPGGLAGAGLGAGLGAFFAVTFPGALVPGGVADAAAAYK AAKAGAGLGGVPGVGGGLGVSAGAVVPQPGAGVKPGKVPVGLPGVYPGGVLPGARFPGVGVLPGVPTGAGVKPKAPG VGGAFAFAGIPGVGPFGGPQP

靶点信息

研究背景	This gene encodes a protein that is one of the two components of elastic fibers. The encoded protein is rich in hydrophobic amino acids such as glycine and proline, which form mobile hydrophobic regions bound by crosslinks between lysine residues. Deletions and mutations in this gene are associated with supravalvular aortic stenosis (SVAS) and autosomal dominant cutis laxa. Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	2006
基因名	ELN
Swiss	P15502
别名	ELN;SVAS;WBS;WS;ADCL1;elastin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn