

VEGFC Rabbit pAb

货号: B11956

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens WB: Homo sapiens , Oryctolagus cuniculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	46kDa
实测分子量	47KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,K-562(Negative control)
细胞定位	Secreted
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human VEGFC (NP_005420.1).
序列	TISFANHTSCRCMSKLDVYRQVHSIIRSLPATLPQCQAANKTCPTNYMWNNHICRLAQEDFMFSSDAGDDSTDGFHDI CGPNKELDEETCQCVCAGLR

靶点信息

研究背景	The protein encoded by this gene is a member of the platelet-derived growth factor/vascular endothelial growth factor (PDGF/VEGF) family. The encoded protein promotes angiogenesis and endothelial cell growth, and can also affect the permeability of blood vessels. The proprotein is further cleaved into a fully processed form that can bind and activate VEGFR-2 and VEGFR-3 receptors.
基因ID	7424

基因名	VEGFC
Swiss	P49767
别名	VEGFC;Flt4-L;LMPH1D;VRP

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn