

Adiponectin Rabbit pAb

货号: B11958

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Sus scrofa , Mesocricetus auratus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	26kDa
实测分子量	30KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Rat heart
细胞定位	Secreted
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 19-104 of human Adiponectin (NP_004788.1).
序列	ETTTQGPVLLPLPKGACTGWMAGIPGHPGHNGAPGRDGRDGTGPEKGEKGDPLIGPKGDIGETGVPGAEGPRGFPGIQGRKGEP

靶点信息

研究背景	This gene is expressed in adipose tissue exclusively. It encodes a protein with similarity to collagens X and d VIII and complement factor C1q. The encoded protein circulates in the plasma and is involved with metabolic and hormonal processes. Mutations in this gene are associated with adiponectin deficiency. Multiple alternatively spliced variants, encoding the same protein, have been identified.
基因ID	9370

基因名	ADIPOQ
Swiss	Q15848
别名	ADIPOQ;ACDC;ACRP30;ADIPQTL1;ADPN;APM-1;APM1;GBP28

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn