

EPHX2 Rabbit pAb

货号: B11981

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	55kDa/57kDa/62kDa
实测分子量	63kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,LO2,SW480,Mouse liver,Mouse kidney,Mouse heart,Rat kidney
细胞定位	Cytoplasm,Peroxisome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-240 of human EPX2 (NP_001970.2).
序列	LALPRGLLNDAFQKGGPEGATTRLMKGEITLSQWIPLMEENCRKCSETAKVCLPKNFSIKEIFDKAISARKINRPMLQAALMLRKKGFTTAILTNTWLDDRAERDGLAQLMCELKMHFDLIESCQVGMVKPEPQIYKFLDLTKASPSEVVFLDDIGANLKPARDLGMVTILVQDQDTALKELEKVTGIQLLNTAPLPTSCNPSPDMSHG

靶点信息

研究背景	This gene encodes a member of the epoxide hydrolase family. The protein, found in both the cytosol and peroxisomes, binds to specific epoxides and converts them to the corresponding dihydrodiols. Mutations in this gene have been associated with familial hypercholesterolemia. Alternatively spliced transcript variants have been described.
------	--

基因ID	2053
基因名	EPHX2
Swiss	P34913
别名	EPHX2;CEH;SEH

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn