

HGS Rabbit pAb

货号: B11983

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Rattus norvegicus IB: Homo sapiens
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:500 - 1:1000
理论分子量	76kDa/86kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1,OVCAR3,MCF7,HepG2,Mouse testis,Rat testis
细胞定位	Cytoplasm,Cytoplasmic side,Early endosome membrane,Endosome,Peripheral membrane protein,multive sicular body membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 548-777 of human HGS (NP_004703.1).
序列	EQQKQTVQMRAQMPAFPLPYAQLQAMPAAGGVLYQPSGPASFSTFSPAGSVEGSPMHGVYMSQPAPAAGPYPSMPSTAADPSMVSAYMYPAGATGAQAAPQAQAGPTASPAYSSYQPTPTAGYQNVASQAPQSLPAISQPPQSSTMGYMGSQSVSMGYQPYNMQNLMTTLP SQDASLPPQQPYIAGQQPMYQQMAPSGGPPQQPPVAQQPQAQGPPAQGSEAQLISFD

靶点信息

研究背景	The protein encoded by this gene regulates endosomal sorting and plays a critical role in the recycling and degradation of membrane receptors. The encoded protein sorts monoubiquitinated membrane proteins into the multivesicular body, targeting these proteins for lysosome-dependent degradation.
基因ID	9146
基因名	HGS
Swiss	O14964
别名	HGS;HRS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn