

E2F7 Rabbit pAb

货号: B11999

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	80kDa/98kDa/99kDa
实测分子量	105kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SP2/0,Mouse spleen, Rat spleen
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 540-740 of human E2F7 (NP_976328.2).
序列	LLAGQPLVYVPSASLFMLYGSLQEGPASGSGSERDDRSSEAPATVELSSAPSAQKRLCEERKPQEEDEPATKRQSREYED GPLSLVMPKKPSDSTDLASPKTMGNRASIPKDIHVNGQLPAAEEISGKATANSLVSSEWGNPSRNTDVEKPSKENESTKE PSLLQYLCVQSPAGLNGFNVLLSGSQTPPTVGPSSGQLPS

靶点信息

研究背景	E2F transcription factors, such as E2F7, play an essential role in the regulation of cell cycle progression (Di Stefano et al., 2003 [PubMed 14633988]).
基因ID	144455

基因名	E2F7
Swiss	Q96AV8
别名	E2F7

产品验证



Western blot analysis of E2F7 expressed in Mouse spleen,Rat spleen,HeLa using E2F7 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. B locking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn