

FAM83A Rabbit pAb

货号: B12000

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IP: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	32kDa/40kDa/47kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse pancreas
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human FAM83 A (NP_116288.2).
序列	MSRSRHLGKIRKRLQDVKSQWVRPARADFS DNESARLATDALLDGGSEAYWRVLSQEGEVDFLSSVEAQYIQAQAREPPCPPDTLGGAEAGPKGLDSSSLQSGTYFPVASEGSEPALLSWASAEKPYLKEKSSATVYFQTVKHNNIRD LVRR C ITRTSQ

靶点信息

研究背景	Probable proto-oncogene that functions in the epidermal growth factor receptor/EGFR signaling pathway. Activates both RAS/MAPK and PI3K/AKT/TOR signaling cascades downstream of EGFR. Required for the RAS/MAPK signaling cascade activation upon EGFR stimulation, it also activates both signaling cascades independently of EGFR activation.
基因ID	84985

基因名	FAM83A
Swiss	Q86UY5
别名	FAM83A;BJ-TSA-9

产品验证



Western blot analysis of FAM83A expressed in Mouse lung,Mouse pancreas using FAM83A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. B locking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn