

FSHR Rabbit pAb

货号: B12003

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Gallus gallus IF: Gallus gallus WB: Gallus gallus , Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	71kDa/75kDa/78kDa/81kDa
实测分子量	72KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,Mouse testis,Rat ovary
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 631-695 of human FSHR (NP_000136.2).
序列	KNFRRDFILLKCGCYEMQAQIYRTETSSTVHNTHPRNGHCSSAPRVTSGSTYILVPLSHLAQN

靶点信息

研究背景	The protein encoded by this gene belongs to family 1 of G-protein coupled receptors. It is the receptor for follicle stimulating hormone and functions in gonad development. Mutations in this gene cause ovarian dysgenesis type 1, and also ovarian hyperstimulation syndrome. Alternative splicing results in multiple transcript variants.
------	--

基因ID	2492
基因名	FSHR
Swiss	P23945
别名	FSHR;FSHR1;FSHRO;LGR1;ODG1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn