

PIK3C3/VPS34 Rabbit pAb

货号: B12024

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IP: Nicotiana tabacum WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:100 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	101kDa
实测分子量	100KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,C6,Mouse lung,Mouse brain,Rat brain
细胞定位	Cytoplasmic vesicle,Late endosome,Midbody,autophagosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 130-412 of human PIK3 C3/VPS34 (NP_002638.2).
序列	FGKYGMFRQGMHDLKVWPNVEADGSEPTKTPGRSTSSTLSEDDQMSRLAKLTKAHRQGHMVKVDWLDRLTFREIEMINES EKRSSNFMYLMVEFRCVKCDDKEYGIVYYEKDGEDESSPILTSFELVKVPDPQMSMENLVESKHHKLARSLRSGPSDHDLK PNAATRDQLNIIIVSYPPTKQLTYEEQDLVWKFRYYLTNQEKALTFLKCVNWDLPQEAQALELLGKWKPMDVEDSLELLS SHYTNPTVRRYAVARLRQADDEDLLMYLLQLVQALKYENFDDIK

靶点信息

研究背景	Catalytic subunit of the PI3K complex that mediates formation of phosphatidylinositol 3-phosphate; different complex forms are believed to play a role in multiple membrane trafficking pathways: PI3KC3-C1 is involved in initiation of autophagosomes and PI3KC3-C2 in maturation of autophagosomes and endocytosis. As part of PI3KC3-C1, promotes endoplasmic reticulum membrane curvature formation prior to vesicle budding. Involved in regulation of degradative endocytic trafficking and required for the abscission step in cytokinesis, probably in the context of PI3KC3-C2. Involved in the transport of lysosomal enzyme precursor to lysosomes. Required for transport from early to late endosomes (By similarity).
基因ID	5289
基因名	PIK3C3
Swiss	Q8NEB9
别名	PIK3C3;VPS34;Vps34;hVps34

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn