

FABP7 Rabbit pAb

货号: B12036

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IF: Mus musculus , Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:100 - 1:500 IF/ICC: 1:50 - 1:200
理论分子量	14kDa/18kDa
实测分子量	13KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,A-549
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human FABP7 (NP_001437.1).
序列	MVEAFCATWKLTN SQNFDEYMKALGVGFATRQVGNVTKPTVIISQEGDKVVIRTLSTFKNTEISFQLGEEFDETTADDRNC KSVVSLDGDGLVHIQKWDGKETN FVREIKDGKMVMTLTFGDVVAVRHYEKA

靶点信息

研究背景	The gene encodes a small, highly conserved cytoplasmic protein that bind long-chain fatty acids and other hydrophobic ligands. The encoded protein is important in the establishment of the radial glial fiber in the developing brain. Alternative splicing and promoter usage results in multiple transcript variants encoding different isoforms. Pseudogenes of this gene are found on multiple chromosomes.
------	--

基因ID	2173
基因名	FABP7
Swiss	O15540
别名	FABP7;B-FABP;BLBP;FABPB;MRG

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn