

p47phox/NCF1 Rabbit pAb

货号: B12038

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Rattus norvegicus IHC: Homo sapiens,Mus musculus
应用	WB IF/ICC IP
推荐浓度	WB: 1:100 - 1:500 IF/ICC: 1:50 - 1:100 IP: 1:50 - 1:100
理论分子量	30kDa/44kDa
实测分子量	47KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,Daudi,Jurkat(Negative control)
细胞定位	Cytoplasm,Cytoplasmic side,Membrane,Peripheral membrane protein,cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-390 of human p47phox/NCF1 (NP_000256.4).
序列	MGDTFIRHIALLGFEKRFVPSQHYVYMFLVKWQDLSEKVVYRRFTEIYEFHKTLEMFPIEAGAINPENRIIPHLPAKWFDG QRAAENRQGTLT EYCSLMSLPTKISRCPHLLDFFKVRPDDLKLPDNTQTKKPETYLMPKDGKSTATDITGP IILQTYRAIAN YEKTS GSEMA LSTGDVVEVVEKSESGWWFCQMKAKRGWIPASFLEPLDSPDETEDPEPNYAGEPYVAIKAYTAVEGDEV SLLEGEAVEVIHKLLDGWWVIRKDDVTGYFPSMYLQKSGQDV SQAQRQIKRGAPPRRSSIRNAHSIHQRSRKRLSQDAYR RNSVRFLQQRRRQARPGPQSPGSPLEEEERQTQRSKPQPAVPPRPSADLILNRCSESTKRKLASAV

靶点信息

研究背景	The protein encoded by this gene is a 47 kDa cytosolic subunit of neutrophil NADPH oxidase. This oxidase is a multicomponent enzyme that is activated to produce superoxide anion. Mutations in this gene have been associated with chronic granulomatous disease.
基因ID	653361
基因名	NCF1
Swiss	P14598
别名	NCF1;NCF1A;NOXO2;SH3PXD1A;p47phox

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn