

F4/80 Rabbit pAb

货号: B12122

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Mus musculus , Ovis aries
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	75-350KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RAW264.7,Neuro-2a(Negative control)
细胞定位	external side of plasma membrane,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 390-641 of mouse F4/80 (NP_034260.1).
序列	VQGRDLGYASFCTLVNATFTILDNTCENKSAPVSLQSAATSVSLVLEQATTWFELSKEETSTLGITILLETVESTMLAALLIPS GNASQMIQTEYLDIESKVINEECKENESINLAARGDKMNVGCFIIKESVSTGAPGVAFVSFAHMESVLNERFFEDGQSFRKL RMNSRVVGGTVTGEKKEDFSKPIIYTLQHIQPKQKSERPICVSWNTDVEDGRWTPSGCEIVEASEHTVCSCNRMANLAI MASGEL

靶点信息

研究背景	Orphan receptor involved in cell adhesion and probably in cell-cell interactions specifically involving cells of the immune system. May play a role in regulatory T-cells (Treg development.
基因ID	13733

基因名	
Swiss	Q61549
别名	DD7A5-7; EGF-TM7; Emr1; F4/80; Gpf480; Ly71; TM7LN3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn