

15-PGDH/HPGD Rabbit pAb

货号: B12168

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Human breast cancer , Homo sapiens WB: Mus musculus , Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	15kDa/19kDa/21kDa/28kDa
实测分子量	25KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,THP-1(Negative control)
细胞定位	Cytoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of human 15-PGDH/15-PGDH/HPGD (NP_000851.2).
序列	MHVNGKVALVTGAAQGIGRAFAEALLLKGA KVALVDWNLEAGVQCKAALDEQFEPQKTLFIQCDVADQQQLRDTFRKVV DHFGRLDILVNNAGVNNEKNWEKTLQINLVVISGTYLGLDYMSKQNGGEGGIIINMSSLAGLMPVAQQPVYCASKHGIV GFTRSAALAANLMNSGVRLNAICPGFVNTAILESIEKEENMGQYIEYKDH IKDMIKYYGILDPLIANG LITLIEDDALNGAIMKI TTSKGIHFQDYDTTPFQAKTQ

靶点信息

研究背景	This gene encodes a member of the short-chain nonmetalloenzyme alcohol dehydrogenase protein family . The encoded enzyme is responsible for the metabolism of prostaglandins, which function in a variety of physiologic and cellular processes such as inflammation. Mutations in this gene result in primary autosomal recessive hypertrophic osteoarthropathy and cranioosteoarthropathy. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	3248
基因名	HPGD
Swiss	P15428
别名	HPGD;15-PGDH;PGDH;PGDH1;PHOAR1;SDR36C1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn