

GABRA2 Rabbit pAb

货号: B12178

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: mouse brain , Rattus norvegicus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	51kDa/52kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1,Mouse brain
细胞定位	Cell junction,Cell membrane,Multi-pass membrane protein,postsynaptic cell membrane,synapse
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human GABRA2 (NP_000798.2).
序列	ARNSLPKVAYATAMDWFIACVAFVFSALIEFATVNYFTKRGWAWD GKSVVNDKKKEKASVMIQNNAYAVAVANYAPNL SKDPVLSTISKSATTPEPNKKP

靶点信息

研究背景	GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA-A receptors, which are ligand-gated chloride channels. Chloride conductance of these channels can be modulated by agents such as benzodiazepines that bind to the GABA-A receptor. At least 16 distinct subunits of GABA-A receptors have been identified. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	2555
基因名	GABRA2
Swiss	P47869
别名	GABRA2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn