

Phospho-MAPKAPK-2/MK2-T334 Rabbit pAb

货号: B12208

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	IHC
推荐浓度	IHC: 1:50 - 1:100
理论分子量	42kDa/45kDa
实测分子量	42kDa/45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	A synthetic phosphorylated peptide around T334 of human MAPKAPK2 (NP_116584.2?).
序列	PQTPL

靶点信息

研究背景	This gene encodes a member of the Ser/Thr protein kinase family. This kinase is regulated through direct phosphorylation by p38 MAP kinase. In conjunction with p38 MAP kinase, this kinase is known to be involved in many cellular processes including stress and inflammatory responses, nuclear export, gene expression regulation and cell proliferation. Heat shock protein HSP27 was shown to be one of the substrates of this kinase in vivo. Two transcript variants encoding two different isoforms have been found for this gene.
基因ID	9261

基因名	MAPKAPK2
Swiss	P49137
别名	MAPKAPK2;MAPKAP-K2;MK-2;MK2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn