

ECSIT Rabbit pAb

货号: B12228

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Mus musculus WB: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	24kDa/33kDa/49kDa
实测分子量	49kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	Cytoplasm,Mitochondrion,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-296 of human ECSIT (NP_001135936.1).
序列	SEQSLVPSPEPRQRPTKALVPFEDLFGQAPGGERDKASFLQTVQKFAEHSVRKRGHIDFIYLALRKMREYGVERDLAVYN QLLNIFPKEVFRPRNIIQRIFVHYPRQCEGIAVLEQMENHGVMPNKETEFLLIQIFGRKSYPMCLKLVRCLKWFPRFMNVNPFPPVRDLPQDPVELAMFGLRHMEPDLSARVTIYQVPLPKDSTGAADPPQPHIVGSGRDAGGVEPLLPDAAGPGVCEEWLGG LRV

靶点信息

研究背景	Adapter protein of the Toll-like and IL-1 receptor signaling pathway that is involved in the activation of NF -kappa-B via MAP3K1. Promotes proteolytic activation of MAP3K1. Involved in the BMP signaling pathway. Required for normal embryonic development (By similarity).
------	---

基因ID	51295
基因名	ECSIT
Swiss	Q9BQ95
别名	ECSIT;SITPEC

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn