

TSSK4 Rabbit pAb

货号: B12230

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	Co-IP: Mus musculus WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa/37kDa/38kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis
细胞定位	sperm flagellum
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 79-328 of human TSSK4 (NP_777604.2).
序列	VLRHKYLINFYRAIESTSRVYIILELAQGGDVLEWQRYGACSEPLAGKWFSQTLGIAYLHKSIVHRDLKENLLLDKWENV KISDFGFAKMVPSNQPVGCSPSYRQVNCFSHLSQTYCGSFAYACPEILRGLPYNPFLSDTWSMGVILYTLVVAHLPFDDTN LKLLRETQKEVTFPANHTISQECKNLILQMLRQATKRATILDIIKDSWVLKFQPEQPTHEIRLLEAMCQLHNTTKQHQSQTIT T

靶点信息

研究背景	This gene encodes a member of the testis-specific serine/threonine kinase family. The encoded protein is thought to be involved in spermatogenesis via stimulation of the CREB/CRE responsive pathway through phosphorylation of the cAMP responsive element binding protein transcription factor. Alternative splicing results in multiple transcript variants.
------	--

基因ID	283629
基因名	TSSK4
Swiss	Q6SA08
别名	TSSK4;C14orf20;STK22E;TSK-4;TSK4;TSSK5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn